




Default Discover:Board Templates and Setups

Discover:Board Templates (Sorted Alphabetically)

Template Name	Picture	Description
2 Keys.tmb		Blank Discover:Board Template
4 Keys.tmb		Blank Discover:Board Template
4 Separate Keys.tmb		Blank Discover:Board Template
8 Keys 2 Row X 4 Col.tmb		Blank Discover:Board Template






Default Discover:Board Templates and Setups

Madentec Limited, 4664 – 99th Street, Edmonton, Alberta, Canada T6E 5H5

(877) 623-3682 • www.madentec.com

Enabling Extraordinary People to do Ordinary Things!



8 Separate Keys.tmb		Blank Discover:Board Template
15 Separate Keys.tmb		Blank Discover:Board Template
32 Keyguard Template.tmb		Template for creating overlays that fit the 32 hole keyguard. Same exact template as the 32 Keys 4 Row X 8 Col.tmb – just renamed for ease of selection.
32 Keys 4 Row X 8 Col.tmb		Blank Discover:Board Template
32 Separate Keys.tmb		Blank Discover:Board Template






Default Discover:Board Templates and Setups

Madentec Limited, 4664 – 99th Street, Edmonton, Alberta, Canada T6E 5H5

(877) 623-3682 • www.madentec.com

Enabling Extraordinary People to do Ordinary Things!



64 Keyguard Template.tmb		<p>Template for creating overlays that fit the 64 hole keyguard. Same exact template as the 64 Keys 8 Col X 8 Row.tmb – just renamed for ease of selection.</p>
64 Keys 8 Col X 8 Row.tmb		<p>Blank Discover:Board Template</p>
128 Keyguard Template.tmb		<p>Template for creating overlays that fit the 128 hole keyguard. Same exact template as the 128 Keys 16 Col X 8 Row.tmb – just renamed for ease of selection.</p>
128 Keys 16 Col X 8 Row.tmb		<p>Blank Discover:Board Template</p>
ABC.tmb		<p>This is a mostly completed template, providing the ability to quickly make a variety of setups. All of the routine keys of the keyboard, as well as mouse emulation keys, are completed to save a great deal of editing time. There are several blank keys to customize to whatever program/purpose. An unlimited number of setups may be created based on this template.</p>




Default Discover:Board Templates and Setups

Madentec Limited, 4664 – 99th Street, Edmonton, Alberta, Canada T6E 5H5

(877) 623-3682 • www.madentec.com

Enabling Extraordinary People to do Ordinary Things!



Activity 1.tmb		Blank Discover:Board Template
Activity 2.tmb		Blank Discover:Board Template
Activity 3.tmb		Blank Discover:Board Template
Activity 4.tmb		Blank Discover:Board Template
Activity 5.tmb		Blank Discover:Board Template



Default Discover:Board Templates and Setups

Madentec Limited, 4664 – 99th Street, Edmonton, Alberta, Canada T6E 5H5

(877) 623-3682 • www.madentec.com

Enabling Extraordinary People to do Ordinary Things!



Activity 6.tmb		Blank Discover:Board Template
Qwerty.tmb		<p>This is a mostly completed template, providing the ability to quickly make a variety of setups. All of the routine keys of the keyboard, as well as mouse emulation keys, are completed to save a great deal of editing time. There are several blank keys to customize to whatever program/purpose. An unlimited number of setups may be created based on this template.</p>

Default Discover:Board Templates and Setups






Madentec Limited, 4664 – 99th Street, Edmonton, Alberta, Canada T6E 5H5

(877) 623-3682 • www.madentec.com

Enabling Extraordinary People to do Ordinary Things!



Discover:Board Setups (Sorted Alphabetically)

Setup Name	Picture	Description
Advanced Qwerty.sub		Expanded keyboard layout with complete mouse controls
Easy ABC.sub		Basic word processing and numbers programs; contains upper-case letters on the setup
Easy Mouse.sub		Any application that requires simple mouse movements
Easy Numbers.sub		Basic numbers programs
Easy Qwerty.sub		Standard keyboard layout with complete mouse controls






Default Discover:Board Templates and Setups

Madentec Limited, 4664 – 99th Street, Edmonton, Alberta, Canada T6E 5H5

(877) 623-3682 • www.madentec.com

Enabling Extraordinary People to do Ordinary Things!



Internet Explorer.sub		Internet Explorer access
Mouse Complex Large.sub		A more complex mouse overlay that is large for those who do not have very good fine motor control
Mouse Complex Small.sub		A more complex mouse overlay that is small for those who have a smaller range of motion
Mouse Simple Large.sub		A simple mouse overlay that is large for those who do not have very good fine motor control
Mouse Simple Small.sub		A simple mouse overlay that is small for those who have a smaller range of motion




Default Discover:Board Templates and Setups

Madentec Limited, 4664 – 99th Street, Edmonton, Alberta, Canada T6E 5H5

(877) 623-3682 • www.madentec.com

Enabling Extraordinary People to do Ordinary Things!



<p>Mouse.sub</p>		<p>A simple mouse overlay that is large for those who do not have very good fine motor control and includes mouse clicking functions</p>
<p>Outlook Express.sub</p>		<p>Outlook Express access</p>
<p>Yes No.sub</p>		<p>A simple setup used in evaluating if a student can answer yes/no questions. Easily adapted for choice making activities, etc.</p>

Default Discover:Board Templates and Setups

Madentec Limited, 4664 – 99th Street, Edmonton, Alberta, Canada T6E 5H5

(877) 623-3682 • www.madentec.com

Enabling Extraordinary People to do Ordinary Things!

